

Fibertec™ 1023 Dietary Fibre Solutions



Tecator™ Line

The Fibertec™ 1023 Dietary Fibre Analysis System is a semi-automated incubation and filtration system for quantitative determination of dietary fibre in processed food, cereal products, fruits and vegetables, health food, plant materials etc., using established enzymatic methods, for use in the laboratory.

| Sample | Parameters |
|---------------------------|--|
| Food and Food Ingredients | Total, soluble and insoluble dietary fibre |



Rapid Dietary Fibre analysis

The Fibertec™ 1023 is designed for rapid and rational determination of dietary fibre according to established enzymatic methods, including those approved by AOAC, AACC, NMKL and Asp.

The system includes a shaking water bath and a filtration module for quantitative determination of dietary fibre in a variety of sample types. The filtration module filters and collects six sample solutions and includes a system for rapid dehydration. The filtration step is speeded up through a built-in "Pressure Mode", to break up clogged filter residues during filtration. Filtration time varies from 2-4 minutes for 100 ml of digested solution. The WB 1024 Water bath incubates 12 samples in each batch. The temperature range is from ambient up to 80°C and it is possible to preset three temperatures.

Batch handling processes up to 6 samples at a time saving time. Laboratories around the globe have invested in Fibertec™ instruments that go on working year after year, making Fibertec™ a solution that will keep your lab up to speed, well into the future.

Effective laboratory operations

The Fibertec™ 1023 takes fibre analysis to new levels of efficiency with facilities to match today's demands. The Fibertec™ systems handle samples in specially designed incubation flasks, along with a range of innovative features to make your laboratory operations even simpler:

- Free choice of methods like AOAC, AACC, NMKL and Asp
- Determines total, soluble and insoluble dietary fibre

- Rapid filtration of six samples simultaneously
- Specially designed incubation flasks
- High incubation and separation capacity
- Repeatable readings, S.D. < 0.5% absolute
- Easy and convenient to use
- Saves time, space and manual work
- Minimal downtime and service with long-life components

Laboratory staff will be grateful for the reliability, ease of use and time-saving facilities that make everyday operations simpler and more convenient.

APPLICATIONS

- Total Dietary Fibre (TDF) according to AOAC 991.43 and AACC 32-07
- Soluble (SDF) and Insoluble (IDF) Dietary Fibre according to AOAC 991.42/43 and AACC 32-07
- Total Dietary Fibre (TDF) according to AOAC 985.29 and AACC 32-05
- Soluble (SDF) and Insoluble (IDF) Dietary Fibre according to AOAC 985.29 and AACC 32-05
- Integrated Total Dietary Fibre (ITDF) Codex Definition method according to AOAC 2009.01*

A Total Dietary Fibre Enzyme Kit for 200 determinations (100 duplicate estimations) is available for above. Part Number 60046527

**Note. The Fibertec™ system forms one part of the instrumentation required for the Integrated TDF method.*



Semi-automated incubation and filtration system for quantitative determination of dietary fibre* in food, cereals, fruit and vegetables.
*Note. A Boiling Water Bath, not supplied, is required for estimation of dietary fibre

Technology

The Fibertec™ 1023 determines total dietary fibre as well as the insoluble and soluble fibres separately. It is also useful in sample preparation for detailed studies of fibre fractions, e.g. by chromatographic or colorimetric techniques.

The basic Fibertec™ 1023 includes a WB 1024 Water bath and a Filtration Module*. The WB 1024 Water bath incubates 12 samples in each batch. The temperature range is from ambient to 80°C and it is possible to preset three temperatures. The shaking frequency and amplitude are continuously adjustable.

*A Boiling Water Bath, not supplied, is required for estimation of dietary fibre

The Filtration Module filters and collects 6 sample solutions batch wise and includes a system for rapid dehydration. The filtration step is speeded up through built-in "Pressure Mode". Reversed pressure is applied to facilitate filtration by breaking up clogged filter residues during filtration. Filtration time varies according to sample type, but on average it is 2-4 minutes for 100 ml of digested solution.

The reproducibility (standard deviation) of the Fibertec™ 1023 in determining total, insoluble and soluble dietary fibre is < 0.5% absolute.



Secure your investment with a FossCare™ Support Agreement

Let FOSS take care of you for a maximum return on your analytical investment. Get a four year warranty as part of the new FossCare Premium Preventive Maintenance Agreement or two years as part of any other FossCare agreement. In addition to the peace of mind afforded by the warranty period, the continual preventive maintenance pays off by keeping your analytical instruments working perfectly every day, year after year.

Why preventive maintenance?

As with any analytical solution, it is essential that your FOSS instrument receives regular maintenance to ensure optimal performance and extended lifetime. Avoiding expensive downtime is a matter of following factory standards and preventively replacing parts before they wear out. In turn, this helps ensure reliable and consistent results at the highest level.

Preventive and predictive maintenance combined with global support from 300 dedicated service, application, software and calibration specialists keeps your instrument running perfectly all year round.



Benefits of a FossCare™ Support Agreement:

- Extended Warranty (two or four years depending on the chosen agreement)
- Regular maintenance; the instrument is diagnosed, cleaned, adjusted, tested, fine tuned and recalibrated
- Minimal downtime from replacing components before they are worn out
- Consistent, accurate and reliable results you can always trust
- Preventative maintenance visits when it suits you (your business)
- 24/7 phone support - no need to worry about closing hours or PO
- Low, fixed service budget prevents unexpected expenses
- Discounts on additional services, spares, training, reagents, consumables and software upgrades

Contact your local Foss office for more information.

Technical Specifications

System description

Fibertec™ 1023

- Filtration Module complete, comprising: Filtration Module, Incubation Flasks (4 sets of 6), Standard Crucibles (P2, 2 sets of 6), Basket for Incubation Flasks, 2 Stands for Crucibles, 2 Spray bottles, 2 Water aspiration pumps, Filter aid, Tubing, WB 1024 Water bath, complete, comprising: a tray for 12 Incubation Flasks and a hinged inclined lid.
- Filtration module as above but without WB 1024 Water bath

In addition, we recommend:

- Incubation Flask (set of 6), at least 4 sets
- Crucible, P2 (40 - 100 µm) set of 6, at least 2 sets
- Basket for 12 Incubation Flasks, at least 1
- Total Dietary Fibre Kit for 200 determinations
Note – A Balance, Boiling water bath, Drying oven, Muffle furnace and Kjeldahl system are also required

Accessories:

Alternative filter crucibles:

- 1000 1174 Crucible, P0 (160 - 250 µm), set of 6
- 1000 1171 Crucible, P1 (90 - 150 µm), set of 6
- 1000 1173 Crucible, P3 (16 - 40 µm), set of 6
- 1000 1079 Stand for 6 Crucibles
- CT 193 Cyclotec™ Sample Mill
- CM 190 Cemotec™ Sample Mill
- KN 195 Knifetec™ Sample Mill for high fat, high fibre samples
For protein determination we recommend our Kjelttec analyzers.


Alternative flask trays for other applications, e.g. biochemical work:

- 1000 2203 Tray for E**-flasks, 50 ml
- 1000 2204 Tray for E**-flasks, 100 ml
- 1000 2205 Tray for E**-flasks, 250 ml
- 1000 2206 Tray for E**-flasks, 500 ml
- ** Erlenmeyer

| Performance data: | |
|--------------------------------------|---|
| Filtration Module: | |
| Filtration and dehydration capacity: | 6 samples/batch |
| Filtration rate: | 2-4 minutes/100 ml of digested solution, depending on sample type |
| Repeatability: | S.D. < 0.5% absolute |
| Flask volume: | 600 ml |
| Accuracy: | According to official methods |
| WB 1024 Water bath: | |
| Capacity: | 12 flasks/batch |
| Temperature range: | ambient to 80°C |
| Temperature accuracy: | ± 0.1°C |
| Shaking frequency: | 0-150 strokes/minute |
| Stroke length: | 0-50 mm |

| Installation requirements: | | | | |
|----------------------------|---------------------|-------------------|----------------------|--------|
| Equipment | Power supply | Power consumption | Dimensions w x d x h | Weight |
| 1023 Filtration Module | 210-230 V, 50-60 Hz | 8 | 770 x 280 x 500 | 24 kg |
| WB 1024 Water bath | 210-230 V, 50-60 Hz | 1.500 | 890 x 370 x 320 | 24 kg |

A complete analytical scheme for dietary fibre determination consists of:



FOSS

FOSS
Foss Allé 1
DK-3400 Hilleroed
Denmark

Tel.: +45 7010 3370
Fax: +45 7010 3371

info@foss.dk
www.foss.dk

