

Thermo Scientific Forma 900 Series -86°C Freezers

Daily sample protection and dependability



Rugged construction for everyday use:

- Heavy gauge, cold-rolled steel cabinets with a powder coat finish for a uniform exterior that resists chipping and rust
- 5" (12.7 cm) foamed-in-place, polyurethane insulation
- State-of-the-art refrigeration system improves temperature control, and increases reserve capacity
- Two 1" (2.5 cm) access ports enable use of inexpedient probes
- Vacuum relief port permits easy access after door openings
- Easy-to-remove washable filter provides protection from dust on the condenser, increasing refrigeration performance
- Single-hand operation with padlock-compatible, ergonomic door handle and integrated key lock

Secure inner doors:

- Four inner doors reduce cold air loss and improve temperature recovery after door openings

Centralized information center:

- Microprocessor control and monitoring system ensure that all controls and displays are easy to reach and read
- Power management system with low voltage surge protection and buck/boost
- Optional choice of CO₂ or LN₂ safety back-up systems for additional sample protection in the event of a power or mechanical failure

Thermo Scientific Forma 900 Series Upright Freezers (Temperature Range: -50°C to -86°C) **CFDA Registered**

Catalogue Number	Model Number	Interior	Capacity		Certification	Electrical	Amps/Breaker (Plug) (Plug is country dependent. See chart below)	Max. Shelf Weight	Interior Dimensions	Exterior Dimensions	Shipping Weight
			liters (cu. ft.)	2" cryobox					HxWxD cm (in)	HxWxD cm (in)	kg. (lbs.)
902-ULTS	ULTS1368	Painted Interior	368 (13.0)	240	CFDA / CE	230V / 50Hz	12/15	57 (125)	130.8 x 58.7 x 49.3 (51.5 x 23.1 x 19.4)	197.9 x 91.2 x 83.6 (77.9 x 35.9 x 32.9)	323 (712)
905-ULTS	ULTS1490	Painted Interior	490 (17.3)	320	CFDA / CE	230V / 50Hz	12/15	57 (125)	130.8 x 58.7 x 64.5 (51.5 x 23.1 x 25.4)	197.9 x 91.2 x 98.8 (77.9 x 35.9 x 38.9)	376 (829)
906-ULTS	ULTS1651	Painted Interior	651 (23.0)	400	CFDA / CE	230V / 50Hz	12/15	57 (125)	130.8 x 77.7 x 64.5 (51.5 x 30.6 x 25.4)	197.9 x 110.2 x 98.8 (77.9 x 43.4 x 38.9)	399 (880)







Options (Field-installed must be installed by a qualified professional)

Description		Model Number
LN₂ Backup System Maintains temperature down to -80°C with liquid nitrogen	Factory Installed	1950447-ULTS
	Field Installed	1950448
CO₂ Backup Maintains temperature down to -70°C with CO ₂	Factory Installed	1950445-ULTS
	Field Installed	1950446
Ink Chart Recorder 15.2 cm (6 in.) seven day ink recorder, to -150°C. Includes one red felt top pen and box of chart paper (qty. 50)	Factory Installed	201241-ULTS
	Field Installed	201240-ULTS
Inkless Chart Recorder 15.2 cm (6 in.) seven day ink recorder, to -150°C. Includes box of chart paper (qty. 50)	Factory Installed	201249-ULTS
	Field Installed	201248-ULTS
PT100 Sensor for connection to building automation systems	Factory Installed	1517487
	Field Installed	290211

Accessories

Description	Model Number
Stainless Steel Shelf Kit 368L freezer, 4 door	1519002G01
Stainless Steel Shelf Kit 490L freezer, 4 door	1519002G02
Stainless Steel Shelf Kit 651L freezer, 4 door	1519002G03

Plug Configurations (Choose our standard or speciality plugs to configure your freezer to meet your local requirements)

Model No.	Plug Description	Male End
Standard Plugs		
US125V20A	NEMA 5-20	
US230V15A	NEMA 6-15	
EU230V16A	Continental EU	
Specialty Plugs		
IN230V16A	India	
AU230V16A	Australia	
CH230V16A	China	

Thermo Scientific Forma 900 Series Upright Freezers (Temperature Range: -50°C to -86°C) General Purpose

Catalogue Number	Model Number	Interior	Capacity		Certification	Electrical	Amps/Breaker (Plug) (Plug is country dependent. See chart below)	Max. Shelf Weight	Interior Dimensions	Exterior Dimensions	Shipping Weight
			liters (cu. ft.)	2" cryobox							
902GP-ULTS	902GP	*Painted Interior	368 (13.0)	240	CE	230V / 50Hz	12/15	57 (125)	130.8 x 58.7 x 49.3 (51.5 x 23.1 x 19.4)	197.9 x 91.2 x 83.6 (77.9 x 35.9 x 32.9)	323 (712)
909-ULTS	909				UL / CE	208-230V / 60Hz	12/15 (NEMA 6-15)				
905GP-ULTS	905GP	*Painted Interior	490 (17.3)	320	CE	230V / 50Hz	12/15	57 (125)	130.8 x 58.7 x 64.5 (51.5 x 23.1 x 25.4)	197.9 x 91.2 x 98.8 (77.9 x 35.9 x 38.9)	376 (829)
910-ULTS	910				UL / CE	208-230V / 60Hz	12/15 (NEMA 6-15)				
906GP-ULTS	906GP	*Painted Interior	651 (23.0)	400	CE	230V / 50Hz	12/15	57 (125)	130.8 x 77.7 x 64.5 (51.5 x 30.6 x 25.4)	197.9 x 110.2 x 98.8 (77.9 x 43.4 x 38.9)	399 (880)
957-ULTS	957				UL / CE	208-230V / 60Hz	12/15 (NEMA 6-15)				

*Factory-installed stainless steel interior option available







Options (Field-installed must be installed by a qualified professional)

Description	Model Number	
LN₂ Backup System Maintains temperature down to -80°C with liquid nitrogen	Factory Installed	1950447-ULTS
	Field Installed	1950448
CO₂ Backup Maintains temperature down to -70°C with CO ₂	Factory Installed	1950445-ULTS
	Field Installed	1950446
Ink Chart Recorder 15.2 cm (6 in.) seven day ink recorder, to -150°C. Includes one red felt top pen and box of chart paper (qty. 50)	Factory Installed	201241-ULTS
	Field Installed	201240-ULTS
Inkless Chart Recorder 15.2 cm (6 in.) seven day ink recorder, to -150°C. Includes box of chart paper (qty. 50)	Factory Installed	201249-ULTS
	Field Installed	201248-ULTS
PT100 Sensor for connection to building automation systems	Factory Installed	1517487
	Field Installed	290211
*Stainless Steel Interior	Factory Installed	SS-ULTS

Accessories

Description	Model Number
Stainless Steel Shelf Kit 368L freezer, 4 door	1519002G01
Stainless Steel Shelf Kit 490L freezer, 4 door	1519002G02
Stainless Steel Shelf Kit 651L freezer, 4 door	1519002G03

Plug Configurations (Choose our standard or speciality plugs to configure your freezer to meet your local requirements)

Model No.	Plug Description	Male End
Standard Plugs		
US125V20A	NEMA 5-20	
US230V15A	NEMA 6-15	
EU230V16A	Continental EU	
Specialty Plugs		
IN230V16A	India	
AU230V16A	Australia	
CH230V16A	China	



thermoscientific.com/cold

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Australia +61 39757 4300

Austria +43 1 801 40 0

Belgium +32 53 73 42 41

China +800 810 5118 or

+400 650 5118

France +33 2 2803 2180

Germany national toll free 0800 1 536 376

Germany international +49 6184 90 6000

BRCSFORMA900-CFDA 0116

India toll free 1800 22 8374

India +91 22 6716 2200

Italy +39 02 95059 552

Japan +81 3 5826 1616

Netherlands +31 76 579 55 55

New Zealand +64 9 980 6700

Nordic/Baltic/CIS countries

+358 10 329 2200

Russia +7 812 703 42 15

Spain/Portugal +34 93 223 09 18

Switzerland +41 44 454 12 12

UK/Ireland +44 870 609 9203

USA/Canada +1 866 984 3766

Other Asian countries +852 2885 4613

Countries not listed +49 6184 90 6000

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand