

ELECTRONIC BALANCE

Features:

- Touch key Internal Calibration
- ISO Certified
- Easy Setting
- Expanded Piece Counting Function
- Comparator Function
- Formulation Mode
- Power Saving Function
- Windows Direct Communication Function
- Power Lock
- Should have superior technology like Unibloc for superior accuracy, response and durability

Technical Specification:

- ~Capacity : 220g
- ~Minimum Display : 0.1mg
- ~Repeatability (Standard deviation) : ≤ 0.1 mg
- ~Linearity : ± 0.2 mg
- ~Stabilisation Time : Approx 3 Seconds
- ~Operating Temperature and Humidity Limits :5-40°C 20-85%*2
- ~Temperature Coefficient for
- ~Sensitivity (10-30°C) : ± 2 ppm/ °C
- ~Pan Size (mm) approx : $\varnothing 91$
- ~Main Body Dimensions (mm) approx : 213(W) 356(D) 338(H)
- ~Main Body Weight (kg) approx : 6.2
- ~Power Requirement : 12V, 1A
- ~Internal Calibration : Yes