

## **ELECTRONIC BALANCE**

### **Features:**

- ISO Certified
- Very quick response
- Compact body
- High stability
- Various unit conversions
- Analog bar display
- Square pan Interface capabilities
- Piece counting function

### **Technical Specification:**

- ~Capacity: 220g
- ~Minimum display: 0.001g
- ~Repeatability (standard deviation):  $\sigma \leq 0.001g$
- ~Linearity:  $\pm 0.02g$
- ~Response time (approx.): 1.0–1.2 sec
- ~Ambient temperature: 5–40°C / 41–104°F
- ~Temperature coefficient of
- ~Sensitivity (10–35°C):  $\pm 3ppm/^{\circ}C$
- ~Pan size (mm): 100×100 (mm)
- ~Dimensions (approx.) : 170W×240D×75H mm
- ~Weight (approx.): 2.2 kg
- ~Power requirements: AC adapter; 90V–264V, -15%–+10%, 50/60Hz
- ~Standard accessories: Simple windbreak (brown smoke color), Protective in-use cover
- ~Display: LCD