## **Refrigerated Micro Centrifuge**

Max RCF: 21,000 xg or Better Max Speed: 14,000 rpm or Better

Max Noise Level: 56 dBA (ventilated); 50 dBA (refrigerated)

Acceleration/Deceleration Rates: 11/12 seconds (ventilated); 10/12

seconds (refrigerated)

Time Set Range: 1 min - 99 min; 1 min increments + HOLD mode Temperature Range: Set from -9 °C to +40 °C per 1 °C increment Certifications: CSA certified, CE marked, IVD compliant, Certified biocontainment with Biocontainment certification by a Third Party Technical Standards: IEC 61010-1, IEC 61010-2, IEC 61010-020

## Should be supplied with:

24 x 1.5/2.0 mL rotor

with ClickSeal biocontainment lid centified by a 3rd Party

10 x 5 mL rotor

with ClickSeal biocontainment lid certified by a 3rd Party for Nucleic Acid

Centrifugation

Personal Centrifuge to be included